

Equipment

Regulation: #3 RCDSO

"Motors...such components should be covered with barriers that are changed after each patient"

However, a lot of dental offices have reported that it is difficult to cover the motors with barriers, therefore additional motors may be needed.

Air Motor



(#13858)

Regulation: #3 RCDSO

"The handpiece should then be removed and, following cleaning and disinfection of clinical contact surfaces, another sterilized handpiece may be attached for used with the next patient."

Piezo Scaler Handpiece



(#14568)

Regulation: #1 IPAC

"A sharp container must be easily accessible in every "point of use" area (e.g., individual examining room)."

Sharps Container



(#17039)

Regulation: #3 RCDSO

"Critical and semi-critical instruments should be processed in a manner that will maintain sterility during storage. Suitable packaging materials include wrapped perforated instrument cassettes."

Sterilization Cassette



(#14744)

Supplies

Regulation: # 1 IPAC

"Suction tubing shall be disposable."

Saliva Ejectors with Affixed Tips (Latex Free)



(#18879)



(#18880)

Regulation: #1 IPAC

"Equipment and surfaces in the clinical office setting must be cleaned with approved hospital-grade cleaners and disinfectants. Equipment cleaning/disinfection should be done as soon as possible after items have been used."

Ethanol Disinfectant Cleaner



(#11015)

Isopropyl Alcohol 70%



(#10978)

Regulation: # 3 RCDSO

"Suction lines should be purged between patients by aspirating water or an appropriate cleaning solution, suction lines should be cleaned at least once per week."

Evacuation System Washer



(#11008)

Saliva Ejector Suction



(#11682)

Regulation: #1 IPAC

"Clean instrument with an instrument detergent/enzymatic diluted according to the detergent manufacturer's directions or, alternatively, in an ultrasonic washer."

Ultrasonic Tablets
(64Pcs/Box)



(#12035)

Ultrasonic General Purpose Cleaner



(#11005)

HVE Suction



(#11679)



(#11680)

Regulation: #1 IPAC

"For the instruments (semi-critical medical equipment) that contact mucous membranes that cannot tolerate sterilization, High Level Disinfectants must be used."

Activated 2% Glutaraldehyde



(#10976)

Legend

*IPAC: Public Health Ontario Guidelines *RCDSO: Royal College of Dental Surgeon of Ontario Guidelines (Please refer to the guideline list we have sent earlier.)