

STERILIZATION QUICK GUIDE

ENSURE PROPER USE OF YOUR STERILIZER. FOLLOW THESE STEPS FOR STERILITY ASSURANCE, **EVERY TIME**

STEAM STERILIZERS

Steam sterilizers are recommended as the sterilization method of choice by the CDC.

Gravity Displacement & SFPP Sterilizer

MIDMARK: M3 / M9 / M11

SCICAN: Statim 2000 / 5000 / 2000 G4 / 5000 G4

TUTTNAUER™: EZ9™, EZ10™, EZ10k™, 3870EA, EZ9 Plus, EZ11 Plus, 2340M, 2540M, 2540MK, 3870M, Valueklave™ 1730

1.



Class 4 Pouches

For use with **EVERY LOAD**
(Internal/external multi-parameter
chemical indicators)

2.



Steam™ Type 5* Integrator

For use with **EVERY PACK** or at least
EVERY LOAD for added
sterility assurance

3.



ConFirm® Biological Monitors

Spore test **AT LEAST DAILY** for each
sterilizer per cycle type

Dynamic Air Removal (Class B) Sterilizer

Class B – Requires daily air removal test (Bowie-Dick) to verify vacuum pump is working properly.

SCICAN: Bravo

W & H: Lisa

TUTTNAUER: Elara 11

1.



Run Warm-up Cycle

with an empty chamber at
the start of **EVERY DAY**

2.



AirView™ II Bowie-Dick Test

Use in empty chamber after
warm-up cycle **EVERY DAY**

3.



Sure-Check® Pouches

For use with **EVERY LOAD**
(Internal/external multi-parameter
chemical indicators)

4.



Steam™ Type 5* Integrator

For use with **EVERY PACK** or
at least **EVERY LOAD**
for added sterility assurance

5.



ConFirm® Biological Monitors

Spore test **AT LEAST DAILY** for each
sterilizer per cycle type

Order Info.
at Next Page



STEAM STERILIZERS RECOMMENDATION

Sure-Check™ Self-Seal Sterilization Pouches

500/box (11801) 3.5" x 9" Reg. \$35.95

200/box (11805) 2.25" x 4" \$14.89

(11803) 2.75" x 9" \$19.99

(11804) 3.5" x 5.25" \$19.99

(11800) 3.5" x 9" \$19.99

(11802) 5.25" x 6.5" \$19.99

(11799) 5.25" x 10" \$29.95

(11797) 5.25" x 15" \$49.95

(11796) 7.50" x 13" \$49.25

(11798) 8" x 16" \$69.95

100/box (11794) 10" x 15" \$58.37

(11795) 12" x 18" \$64.89



Built-in, multi-parameter, internal/external process indicators for all three sterilization criteria of time, temperature and steam.

AirView™ II Bowie-Dick Test Pack

Easy-to-interpret indicator; smaller size reduces waste; does not contain lead.

30/case Reg. Sale (19169) \$241.99 \$199.99



Steam Plus™ Type 5 Integrators

Immediate-readout; confirm sterilization conditions with a definitive pass or fail result.

(11806) 4" Strips, 100/bag Reg. \$69.99 (19176) Dry Heat Chemical Indicator Strip, 100/box Reg. \$25.15 Sale: \$22.15



ConForm® 10, Result in 10hrs In-Office Biological Monitoring

Easy to use: Simply process indicator in a normal steam sterilization cycle followed by a 10 hour incubation cycle.

(19168) Confirm 10 Starter kit (Incubator, Record Keeper, Manual) Reg. \$628.95 Sale: \$599

(11793) Incubator Machine Reg. \$334.95 Sale: \$319

(11778) Biological Indicators, 25/box Reg. \$143.80 Sale: \$136.95



ConForm® 24, Result in 24hrs In-Office Biological Monitoring

Easy to use: Simply process indicator in a normal steam sterilization cycle followed by a 24 hour incubation cycle.

(19226) Confirm 24 Starter kit (Incubator, 25 Indicators, Crusher, Record Keeper, Manual) Reg. \$481.95 Sale: \$459

(19103) Incubator Machine Reg. \$334.95 Sale: \$319

(19052) Biological Indicators, 25/box Reg. \$125.95 Sale: \$119.99

(19501) SporView Self-Contained Steam BI 100 /box Reg. \$298 Sale: \$268



All product names are trademarks of Crosstex International Inc., a Cantel Medical Company, its affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions. © 2021 Crosstex International, Inc.



Atop Dental Canada Inc. 9-595 Middlefield Rd, Toronto, Ontario, M1V 3S2 Canada Toll Free: 1-877-703-3562 | Fax: 1-866-680-8748 | Email: order@atopdent.ca | Website: atopdent.ca